


LUNA[®] G3

Auto-CPAP

**MAKE PATIENT
COMPLIANCE
A BREEZE.**




React Health's **Luna G3** Auto-CPAP devices come standard with all the features that make sense for you and your patients. Our standard integrated cellular modem makes capturing patient compliance a breeze.

Our standard patient-forward features include integrated heated tubing, pre-heat humidifier function and Auto Ramp offering added comfort to patient therapy as well as a bedside clock and resupply reminders for the convenience of the patient. You will find these features and more, in a smaller footprint, with an easy-to-read full color screen and intuitive interface.

React Health[®] is an American sleep and respiratory device manufacturer and distributor.



DEVICE SPECIFICATIONS



Device Dimensions (L x W x H)	10.4" x 5.7" x 4.5"
Weight	3.75 lbs.
Modes of Operation	Auto-CPAP (LG3600, LG36002)
SD Card	Standard
Pressure Range	4 to 20cmH2O (in 0.5 cm increments)
Ramp Time	0 to 60 minutes
Sound	<28 dBA when device is at 10cmH2O
Water Capacity	360 ml (to maximum fill line)

REPLACEMENT PARTS

LG3600	Luna G3 APAP Device with Cellular Modem-HT
LG36002	Luna G3 APAP Device with Cellular Module without HT
LG34510	Luna G3 Replacement Water Chamber
LG3HT	Luna G3 Integrated Heated Tubing
UN100	Ultra Noir CPAP Standard Tubing (Non-Heated)
LG3F002	Luna G3 Disposable Ultra Fine Filter
LG38008	Luna G3 Reusable Filter
LG3B	Luna G3 Repl Carry Bag
LPG3	Luna G3 Power Supply Adapter
LPAC	Luna G2/G3 Power Cord
PRTG301	Luna G3 Repl Filter Cover- Black
PRTG302	Luna G3 Repl Filter Cover- Ivory

REACT HEALTH OFFERS MORE WAYS TO CONNECT



Cellular



SD Card



QR Code
Scanning

ReactHealth.com

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.
MI.MK.0204, Rev F

