


# Introducing the Luna Family of PAP Devices from React Health

Luna G3 APAP (LG3600, LG36002)



Luna G3 Bilevel 25A (LG3700) and Bilevel ST 30Vt (LG3800)



## LUNA FEATURES:

### The Luna Family of PAP Devices offer:

- RESlex™ pressure relief on exhalation feature for increased patient comfort and compliance.
- Leak compensation to maintain patient's pressure and breathing synchronization.
- Simple release water chamber for easier removal and cleaning.
- Large capacity water chamber holds up to 350ml.
- Large screens and easy to navigate menus for clinicians and patients.
- Integrated heated humidification.
- Event detection for apneas, hypopneas, snoring and central apneas.
- Bilevel ST model offers Target Volume.
- Heated Tubing standard on all Luna G3 Models.
- **FDA Cleared.**

Luna PAP Devices Make Capturing Compliance Easy.

### 4 separate ways to retrieve patient data:



**Cellular\***



**QR Code Scanning**



**SD Card**



**React Health Connect**

Patient data uploads to *React Health Connect*, React Health's Cloud based platform for viewing patient data.

*\*Additional Fees may apply for Cellular*



REACTHEALTH

# PAP COMPARISON CHART



REACTHEALTH

## MANUFACTURER

React Health

React Health

ResMed

Philips



| BRAND                       | LUNA II                    | LUNA G3                    | AIRSENSE® 11                            | DREAMSTATION® 2               |
|-----------------------------|----------------------------|----------------------------|-----------------------------------------|-------------------------------|
| MODES                       | CPAP, Auto CPAP            | CPAP, Auto CPAP            | AutoSet, CPAP, Elite                    | CPAP, Auto CPAP               |
| Pressure Range              | 4-20                       | 4-20                       | 4-20                                    | 4-20                          |
| Exhalation Relief           | Yes (RESlex)               | Yes (RESlex)               | Yes(EPR)                                | Yes (Flex)                    |
| Remote Setting Changes      | Yes                        | Yes                        | Yes                                     | Yes                           |
| Integrated Humidifier       | Yes                        | Yes                        | Yes                                     | Yes                           |
| Sound Level dBa             | 28dBa                      | 28dBa                      | 27dBa                                   | 27dBa                         |
| Compliance Reporting        | Yes                        | Yes                        | Yes                                     | Yes                           |
| Phone App                   | Yes                        | Yes                        | Yes                                     | Yes                           |
| Connectivity Options        | Cellular, QR Code, SD Card | Cellular, QR Code, SD Card | Cellular, SD Card (Optional), Bluetooth | Cellular , SD Card, Bluetooth |
| Heated Tubing               | Optional                   | Standard                   | Optional<br>(subject to availability)   | Optional                      |
| On Screen Resupply Reminder | No                         | Yes                        | No                                      | No                            |

Dreamstation® is a registered trademark of Koninklijke Philips N.V. (Royal Philips)  
AirSense® is a registered trademark of ResMed family of companies.

MI.MK.0223.REV.E.