

React Health PAP Masks: Offering the right fit, the first time

SIESTA FULL FACE MASK



<1%

REFIT RATE AS EVIDENCED BY
SLEEPGLAD
AGGREGATE DATA



SIESTA NASAL MASK

<1%

REFIT RATE AS EVIDENCED BY
SLEEPGLAD
AGGREGATE DATA

RIO II NASAL PILLOWS MASK



<1%

REFIT RATE AS EVIDENCED BY
SLEEPGLAD
AGGREGATE DATA

*Aggregate independent data from SleepGlad Mask Fitter May
2020- April 2021



REACTHEALTH



Our Family Of CPAP Mask Interfaces



REACTHEALTH

React Health
PAP Masks available at:



*React Health is formerly known as 3B Medical, Inc.
MI.MK.0276.REV.A.

The Siesta

Full Face Mask

- No forehead brace allows for increased line of sight - even with your glasses on
- Self-adjusting ergonomic cushion design allows for superior seal
- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- 360-degree dual-swivel elbow for active sleepers and allows for sleep in any position
- Air vent diffuser for a quiet nights sleep
- Anti-phyxia valve safety feature opens for fresh air in case of machine shut off



The Siesta

Nasal Mask

- Padded support arms minimize irritation and marks around cheekbones and eyes
- Quick release clips allow for simple attachment and removal of mask
- Easy change cushion cup snaps into place allows for ease of removal and cleaning
- 360-degree dual-swivel elbow for active sleepers and allows for any sleep position



The Rio™ II

Nasal Pillows Mask

- A lightweight experience- Just 2.3 ounces: One of the lightest masks on the market today
- Dual wall pillows with trampoline bounce for secure fit
- Split back strap allows for a better fit, while Velcro fasteners create added layer of adjustment
- Unique rotating ball and socket swivel for active sleeper

