

# Rio™ II

## Full Face Mask



**The React Health Rio™ II Full Face Mask is a minimal contact mask that minimizes nose discomfort and doesn't interfere with vision or glasses.**



Rio™ II  
Full Face Mask

- ✦ Under-The-Nose Cushion for enhanced comfort by freeing the nasal bridge, which may reduce feelings of claustrophobia.
- ✦ Circular dispersion ventilation system evenly distributes airflow while minimizing noise- as low as 17 Db(A)
- ✦ Modular main frame compatible with all cushion sizes. (S,M,L)
- ✦ Adjustable headgear that provides a comfortable, secure fit for a wide range of head sizes.
- ✦ Features easy release clips- not magnets - that are common to all React Health masks.

# Rio™ II Full Face Mask



*Equipped with innovative magnet-free clips*



*Circular dispersion ventilation system*



*Designed for effortless disassembly, facilitating thorough and efficient cleaning*



www.reacthealth.com info@reacthealth.com  
Phone: 863-226-6285 Fax: 863-226-6284



## MASK PART NUMBERS

**RFF1000** - Rio II Full Face Mask - All Size Kit

**RFF1001** - Rio II Full Face Mask - SM

**RFF1002** - Rio II Full Face Mask -MD

**RFF1003** - Rio II Full Face Mask - LG

## REPLACEMENT SUPPLIES PART NUMBERS

**RFF2001** - Rio II Full Face Repl Headgear

**RFF3001** - Rio II Full Face Repl Cushion - SM

**RFF3002** - Rio II Full Face Repl Cushion - MD

**RFF3003** -Rio II Full Face Repl Cushion- LG

**RFF4001** - Rio II Full Face Mask w/o Headgear - SM

**RFF4002** - Rio II Full Face Mask w/o Headgear - MD

**RFF4003** - Rio II Full Face Mask w/o Headgear - LG

## ReactHealth.com

© 2024, React Health. All rights reserved.  
Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.  
MI.MK.0392, Rev A