

# Siesta 2

## Full Face Mask

The new and improved Siesta 2 Full Face Mask optimizes air flow and offers updated headgear and frame for increased comfort.



**Siesta 2**  
*Full Face Mask*

- ✦ Circular dispersion vent on frame to reduce airflow into user's eyes and minimize disturbance of bed partner
- ✦ Updated forehead-free frame design for improved line of vision
- ✦ Crown-style headgear with shorter and lower straps for a better fit.
- ✦ Features easy release clips- not magnets - that are common to all React Health masks.

# Siesta 2

## Full Face Mask



*Equipped with innovative magnet-free clips*



*Circular dispersion ventilation system*



*Designed for effortless disassembly, facilitating thorough and efficient cleaning*



### MASK PART NUMBERS

SFF21000: Siesta 2 Full Face Mask - All Size Kit

SFF21001: Siesta 2 Full Face Mask - SM

SFF21002: Siesta 2 Full Face Mask -MD

SFF21003: Siesta 2 Full Face Mask - LG

### REPLACEMENT SUPPLIES PART NUMBERS

S22001: Siesta 2 Full Face Repl Headgear

SFF23001: Siesta 2 Full Face Repl Cushion - SM

SFF23002: Full Face Repl Cushion - MD

SFF23003: Full Face Repl Cushion- LG

SFF24001: Full Face Mask w/o Headgear - SM

SFF24002: Full Face Mask w/o Headgear - MD

SFF24003: Full Face Mask w/o Headgear - LG



www.reacthealth.com info@reacthealth.com  
Phone: 863-226-6285 Fax: 863-226-6284

ReactHealth.com

React Health. All rights reserved. Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change. MI.MK.0397, Rev A