

HomeFill®

Oxygen System



NO DELIVERIES.
NO WORRIES.



The **HomeFill®** Oxygen System by React Health promotes patient ambulation by allowing them to fill their own cylinders from their 5-liter or 10-liter oxygen concentrators, all from the comfort of their own homes.



PROVIDER BENEFITS

- ✦ Significantly reduces operating costs through fewer cylinder deliveries and service calls and lower storage requirements, all without sacrificing patient care
- ✦ No extra equipment needed because patients can breathe off concentrator while filling cylinders
- ✦ Encourages compliance with oxygen therapy as patients have no worries about running out of oxygen

PATIENT BENEFITS

- ✦ Eliminating oxygen deliveries increases independence, convenience and comfort
- ✦ Small, lightweight cylinders allow for greater mobility than larger, traditional oxygen cylinders
- ✦ Using small, refillable cylinders indoors can help avoid trips and falls associated with using a concentrator with long oxygen tubing

FEATURES

- ✦ Provides an unlimited, refillable supply of ambulatory oxygen
- ✦ Patient breathes from a continuous flow oxygen concentrator at home
- ✦ Patient protection ensures oxygen purity in the cylinder is always >90%
- ✦ Easy cylinder connections with the OxyLock™ Coupler

SPECIFICATIONS

Dimensions	15" x 20.5" x 16"
Weight	33 lbs.
Electrical Requirements	120VAC +/- 10% @ 60Hz
Limited Warranty	3 Years / 5,000 hours

COMPATIBLE CONCENTRATORS

Platinum® 5NXG
Platinum® 10 Liter
Perfecto® 2 Series



PART NUMBERS

HomeFill Compressor	IOH200
Ready Rack for Platinum® 10	IOH260
Ready Rack for Perfecto® 2 Series	IOH270
Ready Rack for Platinum® 5NXG	IOH280
M9 Cylinder w/ Integrated Conserver	HF2PC9
M9 Cylinder w/ Post Valve	HF2POST9
D Cylinder w/ Post Valve	HF2POSTD



ReactHealth.com

© 2024, Ventec Life Systems, all rights reserved.
Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.
MKT-00101, Rev B