


HomeFill®


Cylinders



REACTHEALTH



React Health offers a variety of HomeFill Cylinders to fit your client's lifestyle and daily oxygen needs. These cylinders range in size and can last up to 15 hours.



 ventec
LIFE SYSTEMS

3 Customizable Options to Meet Patient Needs

HF2PC9



HF2POST9



HF2POSTD



HomeFill Cylinder M9
with Integrated
Conserver

HomeFill Cylinder M9
with CGA870
Compatible Post Valve

HomeFill Cylinder D with
CGA870 Compatible Post
Valve

Specifications

	HF2PC9	HF2POST9	HF2POSTD
Capacity	255	255	422
Fill Time	140 mins	140 mins	3.5 hours
Height	14"	14"	20.8"
Weight	5.1 lbs	5 lbs	6.3 lbs
Diameter	4.4"	4.4"	4.4"
Duration	140mins cont flow	140mins cont flow	3.5 hours cont flow
Pulse setting 1	15.7 hrs	See Conserver Specifications	
Pulse setting 2	10.1 hrs		
Pulse setting 3	7.5 hrs		
Pulse setting 4	6 hrs		
Pulse setting 5	4.9 hrs		

Cylinder Bags

Cylinder Bag for HF2PC9	HF2PC9BAG
Cylinder Bag for HF2POST9	HF2POST9BAG
Cylinder Bag for HF2POSTD	Options Available on Open Market

*duration based on 20 bpm and average pulse volume per setting