

# React Health

Platinum 10  
Oxygen Concentrator



**The Platinum™ 10 High Flow Concentrator is back.**



Experience the epitome of respiratory care with the React Health Platinum 10 Concentrator.

Boasting highly reliable 10-liter flow, seamless HomeFill compatibility, and a sleek one-piece cabinet for effortless cleaning and maintenance.

**The newest addition to our broad line of products to support all of your Respiratory Care needs.**

## PROVIDER BENEFITS

- ✦ Provides continuous flow up to 10LPM
- ✦ HomeFill compatible for inventory flexibility
- ✦ Patient accessible reset button helps minimize emergency service calls
- ✦ May be eligible for QG billing modifier<sup>1</sup>

## PATIENT BENEFITS

- ✦ User friendly design with easy to see control panel, humidifier bottle and filter access
- ✦ HomeFill compatibility allows patients with higher flow prescriptions to access small, lightweight portable oxygen
- ✦ Lower heat generation and power consumption than other high flow concentrators

## PART NUMBERS

Platinum 10 Oxygen Concentrator	<b>IRC10LX02</b>
Gross Particle Cabinet Filter	<b>1107412</b>
Inlet Filter	<b>1131249</b>
Brass Barbed Connector	<b>1180107</b>



## TYPICAL SPECIFICATIONS

Dimensions	26.5" x 18.5" x 14.5"
Weight	53 lbs.
Continuous Oxygen Flow	2 to 10 L/min
Oxygen Concentration	87-95.6%
Oxygen Outlet Pressure	9 PSI +/- 0.5 PSI
Noise Level	<58 dBA <sup>2</sup>
Power Consumption	585W
HomeFill Compatible	Yes
Warranty	3 Years

### ReactHealth.com

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.  
MKT-00116, Rev A

1- Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program.  
2- Average weighted for 10 L/min flow.