

React Health Platinum[®] 10L Oxygen Concentrator



Unlimit Your Oxygen

Helping our Veterans
breathe easier



Enhancing Veterans' quality of life



React Health Platinum® 10L Oxygen Concentrator

Experience the epitome of respiratory care with the React Health Platinum 10L Concentrator. Boasting a highly reliable 10-liter flow, seamless HomeFill compatibility, and a sleek one-piece cabinet for effortless cleaning and maintenance.

- ✦ Provides continuous flow up to 10LPM
- ✦ User friendly design with easy to see control panel, humidifier bottle and filter access.
- ✦ HomeFill compatibility allows patients with higher flow prescriptions to access small, lightweight portable oxygen.
- ✦ Lower heat generation and power consumption than other high flow concentrators.

TYPICAL SPECIFICATIONS

Platinum 10 Oxygen Concentrator	IRC10LXO2
Dimensions	26.5" x 18.5" x 14.5"
Weight	53 lbs.
Continuous Oxygen Flow	2 to 10 L/min
Oxygen Concentration	87-95.6%
Oxygen Outlet Pressure	9 PSI +/- 0.5 PSI
Noise Level	<58 dBA ¹
Power Consumption	585W
HomeFill Compatible	Yes
Warranty	3 Years

1- Average weighted for 10 L/min flow.



[ReactHealth.com](https://www.ReactHealth.com)

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician. Technical specifications may change.

MKT-00133, Rev A