

# PAP COMPARISON CHART



## MANUFACTURER

React Health



React Health



ResMed



Philips



BRAND	LUNA II	LUNA G3	AIRSENSE 11	DREAMSTATION 2
MODES	CPAP, Auto CPAP	CPAP, Auto CPAP	AutoSet, CPAP , Elite	Auto
Pressure Range	4-20	4-20	4-20	4-20
Exhalation Relief	Yes (RESlex)	Yes (RESlex)	Yes (EPR)	Yes (Flex)
Remote Setting Changes	Yes	Yes	Yes	Yes
Integrated Humidifier	Yes	Yes	Yes	Yes
Sound Level dBa	28dBa	26dBa	27dBa	28dBa
Compliance Reporting	Yes	Yes	Yes	Yes
Phone App	Yes	Yes	Yes	Yes
Connectivity	Cellular, Wi-Fi, QR Code, SD Card	Cellular, Wi-Fi, QR Code, SD Card	Cellular, SD Card	Cellular, Wi-Fi , SD Card
Heated Tubing	Optional	Standard	Optional (subject to availability)	Optional
On Screen Resupply Reminder	No	Yes	No	No