



Lumin Wand

Portable UVC Sanitizer

Safely Cleans & Sanitizes Surfaces in Seconds

The Lumin Wand operates using a form of high energy light that works within a specific, narrow wavelength to reduce bacteria.*

The Lumin Wand utilizes UVC light that cleans surfaces such as bedding, packages, keyboards, toys, smartphones, doorknobs, and any high touch areas!

Lab tested to reduce bacteria* in Seconds

* >99.99% Human Corona virus, >99.99% Influenza H1N1, >99.9% Escherichia coli (e. coli), Staphylococcus aureus, Salmonella enterica and Listeria monocytogenes.



REACTHEALTH




120 MINUTES OF PERFORMANCE


ONE-TOUCH CONTROL


CORDLESS, RECHARGEABLE

LWD200 - Lumin Wand Docking Station

LM5000 - Lumin Wand

Product Specifications

Rated Input:	AC100-240V, 50-60HZ
Input Interface:	DC 12.6V 0.8A
LED Wavelength:	254nm
Size:	13.2" x 2" x 1.75" / 335.5 x 50.8 x 430 mm
Weight:	590g / 1.3 lbs.
Battery Type:	Li-ion 18650
Nominal Voltage:	11.1V
Charging Voltage Limitation:	12.6V

Product Features

- Safety Sensor when device is turned upwards
- Built In proximity sensor
- Long Lasting High powered UVC bulb



REACTHEALTH

www.reacthealth.com
info@reacthealth.com

React Health is formerly known as 3B Medical, Inc.