



REACTHEALTH

New Luna G3 CPAP and Auto-CPAP



Luna G3 CPAP
SKU: LG3500

Luna G3 Auto-CPAP
SKU: LG3600

React Health's all new Luna G3 CPAP and Auto-CPAP devices come standard with all the features that make sense for you and your patients. Our standard integrated cellular modem makes capturing patient compliance a breeze.

Our standard patient-forward features include integrated heated tubing, pre-heat humidifier function and Auto Ramp offering added comfort to patient therapy as well as a bedside clock and resupply reminders for the convenience of the patient. You will find these features and more, in a smaller footprint, with an easy-to-read full color screen and intuitive interface.

Device Specifications

Device Dimensions (L x W x H)	10.4" x 5.7" x 4.5"
Weight	3.75 lbs.
Mode of Operation	CPAP (LG3500), CPAP, Auto-CPAP (LG3600)
SD Card	Standard
Pressure Range	4 to 20cmH2O (in 0.5 cm increments)
Ramp Time	0 to 60 minutes
Sound	<26 dBA when device is at 10cmH2O
Water Capacity	350 ml (to maximum fill line)

Replacement Parts

LG 3500	Luna G3 CPAP Device - HT
LG 3600	Luna G3 APAP Device - HT
LG 34510	Luna G3 Replacement Water Chamber
LG 3HT	Luna G3 Integrated Heated Tubing
UN 100	3B Ultra Noir CPAP Standard Tubing (Non-Heated)
LG31020BX**	Luna G3 Disposable Air Filter, 2 per bag, 3 per box (White)
LG38008BX**	Luna G3 Reusable Filter - Box of 12 Filters (Black)
LG3B	Luna G3 Carry Bag
LPG3	Luna G3 Power Supply Adapter
LPAC	Luna G2/G3 Power Cord
PRTG301	Luna G3 Repl Filter Cover- Black



REACTHEALTH

www.reacthealth.com

info@reacthealth.com

*React Health is formerly known as 3B Medical, Inc.

**Luna G3 has interchangeable filters. You may choose either the LG31020BX (A7038) or LG38008BX (A7039) as a filter for the device.