



# Pulse Oximetry



REACTHEALTH

MI.MK.0260.REV.A.



## Product Specifications

- LED Display
- SpO<sub>2</sub> Measuring Range: 0-100%
- PR Measuring Range: 30-250 bpm
- Accuracy +/- 2% (70%-100%) unspecified (<70%) for SpO<sub>2</sub>
- Power supply 1.5V (AAA size) alkaline batteries
- Continuous run time up to 30 hours
  - Low Battery Indicator
- Compact and lightweight
  - 57 x 31 x 32 mm

## Part Numbers

- PO2BLK Pulse Oximeter (Black)
- PO2BLU Pulse Oximeter (Blue)



REACTHEALTH

[www.reacthealth.com](http://www.reacthealth.com)  
[info@reacthealth.com](mailto:info@reacthealth.com)