



Pulse Oximetry



REACTHEALTH



Product Specifications

- LED Display
- SpO₂ Measuring Range: 0-100%
- PR Measuring Range: 30-250 bpm
- Accuracy +/- 2% (70%-100%) unspecified (<70%) for SpO₂
- Power supply 1.5V (AAA size) alkaline batteries
- Continuous run time up to 30 hours
 - Low Battery Indicator
- Compact and lightweight
 - 57 x 31 x 32 mm

Part Numbers

- PO2BLK Pulse Oximeter (Black)
- PO2BLU Pulse Oximeter (Blue)



REACTHEALTH

www.reacthealth.com
info@reacthealth.com